

# TOOTH MOVEMENT PREDICTABILITY

## Per aligner

	ROTATION	EXPANSION	RETRACTION	BUCCOLINGUAL TORQUE	MESIODISTAL AXIAL MOVEMENT	MESIALIZATION	DISTALIZATION	INTRUSION	EXTRUSION
<b>Upper incisors</b>	3° - 4°	0.3 - 0.4 mm	0.3 - 0.4 mm	1°	1°	0.2 - 0.3 mm	0.2 - 0.3 mm	0.2 mm	0.1 - 0.2 mm
<b>Small sized upper lateral incisors</b>	2° - 3°	0.3 mm	0.3 mm	0.5 °	0.5 °	0.1 - 0.2 mm	0.1 - 0.2 mm	0.1 - 0.15 mm	0.1 - 0.15 mm
<b>Lower incisors</b>	3° - 4°	0.3 - 0.4 mm	0.3 - 0.4 mm	1°	1°	0.2 - 0.3 mm	0.2 - 0.3 mm	0.1 - 0.2 mm	0.1 - 0.15 mm
<b>Canines</b>	3°	0.3 - 0.4 mm	0.3 - 0.4 mm	1°	0.5 ° - 1°	0.2 mm	0.2 mm	0.1 - 0.2 mm	0.1 - 0.15 mm
<b>Premolars</b>	1° - 2°	0.3 - 0.4 mm	0.2 - 0.3 mm	0.5 °	0.5 °	0.1 - 0.2 mm	0.2 mm	0.1 mm	0.1 mm
<b>Molars</b>	1°	0.3 - 0.4 mm	0.2 mm	0.5 °	0.5 °	0.1 - 0.2 mm	0.2 mm	0.1 mm	0.1 mm